



USER MANUAL



WE HAVE YOU COVERED

Our passionate team brings you innovative, quality mobility solutions. Our products are reliable, sleek, and designed with simplicity so that you can rest assured that they work when you need them.

Thanks to our out of the box thinking, we are vigilant to meet your mobility needs with our simple and practical solutions. We recognize the important role that safety plays in consumer choices and are proud to offer a line of products that maintain the safety features found in your vehicle. Backed by third party verification for compliance to applicable Federal Motor Vehicle Safety Standards and Canadian Motor Vehicle Safety Standards, our products are engineered for safety and durability.

Our commitment to quality and safety along with our genuine desire to meet your mobility needs makes for a product we are sure you will be proud to own.





SPEEDY-LIFT

USER GUIDE

To get into your vehicle

The **SPEEDY-LIFT** is a product like no other. Speed and ease are the foundations of its design.

Follow the next few steps:

- After transferring out of your wheelchair, open the sliding door to your vehicle and retrieve the hand held pendant.
- Press and hold the **OUT** button until the **SPEEDY-LIFT** has reached its lowest point outside of the vehicle then release the button.
- Simply roll your wheelchair up to the **SPEEDY-LIFT**, press and hold the **IN** button and watch the **SPEEDY-LIFT** do the rest.
- Once your wheelchair is in the vehicle, release the **IN** button, replace the hand held pendant. You can now close the door and get on your way.



USER GUIDE

To get out of your vehicle

Follow the next few steps:

- Open the sliding door of your vehicle and retrieve the hand held pendant.
- Press and hold the **OUT** button and watch the **SPEEDY-LIFT** do the rest.
- Once your wheelchair is outside of the vehicle, roll it away from the **SPEEDY-LIFT**.
- Press and hold the **IN** button until the **SPEEDY-LIFT** is fully stowed and then release the button.
- Place the hand held pendant inside the vehicle and close the door.

WARNING

DO NOT USE A WHEELCHAIR EXCEEDING 95 LBS.



⚠ WARNING: This product can expose you to chemicals including benzene, which is known to the State of California to cause cancer and birth defects or other reproductive harm. For more information, go to www.P65Warnings.ca.gov.

MANUAL BACKUP SYSTEM

Your new **SPEEDY-LIFT** is equipped with a manual backup system that will allow your lift to be stowed manually. **Please contact your dealer for full diagnostic and needed repairs.**

1. Locate the release handle on the side of the **SPEEDY-LIFT**. ***CAUTION: The SPEEDY-LIFT up / down function is spring loaded.** Place your hand on the top of the unit before lowering the release handle.




2. Press and hold down the release handle and push the **SPEEDY-LIFT** back into the vehicle.



BRACKET DESIGNED FOR A MANUAL FOLDABLE WHEELCHAIR



Remove either 'pin' from the end of the upper cylindrical tube and slide it out from within the bracket. For the lower cylindrical tube, you must remove both 'pins'. Once you have done that, slide the cylindrical tube to the left or right until the one side is free, then simply slide out the other side. Your wheelchair is now ready to be folded.



3 YEAR LIMITED WARRANTY

SPEEDY-LIFT

Adapt Solutions # 1 Ltd. (Adapt Solutions), warrants to the original purchaser of a **SPEEDY-LIFT** that the equipment is free from defects in material and workmanship for a period of three years from date of purchase.

During the first year of the warranty, Adapt Solutions will supply the replacement parts as well as a pre-set monetary amount (determined by Adapt Solutions) for the repair if a defect in material or workmanship is discovered. After the initial year of this warranty, only parts and components are covered. This warranty does not cover labour and other services after the initial year. Freight and other related repair charges will be the responsibility of the original purchaser.

The only remedy for a defect in one of Adapt Solutions products shall be the repair or the replacement, at the discretion of Adapt Solutions, of the defective part or component. If repair or replacement is not commercially practical or cannot be timely made, Adapt Solutions may decide to refund the purchase price of the equipment instead of repairing or replacing the original equipment.

In no event shall Adapt Solutions be responsible for indirect, incidental or consequential damages, whether such damages arise from claims based on contract, warranty, tort (including negligence), strict liability or product liability.

All implied warranties, including any warranty of merchantability or fitness for a particular purpose, are limited in their duration to the length of the warranty stated above for the affected component.

This warranty is to the original purchaser only, and excludes product damage due to installation error, product misuse, product abuse, accidents, physical damage, damage in shipment, modifications not made by Adapt Solutions, or repairs undertaken by anyone other than authorized distributors.

TO OBTAIN WARRANTY PARTS OR REIMBURSEMENT

You must follow these procedures:

1. Obtain warranty authorization by calling your local Adapt Solutions dealer or Adapt Solutions directly at 1-866-641-0419.
2. Return the faulty Adapt Solutions component/equipment, freight prepaid, to the address provided by your Adapt Solutions dealer or Adapt Solutions with proof of purchase indicating the date purchased.

Adapt Solutions will pay for shipping back to the purchaser within the continental United States and Canada if a defect in materiel or workmanship is discovered. Return freight and repair charges will be the responsibility of the purchaser if the problem is not covered by warranty.

This warranty gives you specific legal rights, and you may also have other rights that vary from state to state/ province to province. Adapt Solutions specifically does not authorize any person to extend the time or scope of this warranty.

For further information regarding this limited warranty, please contact us by calling 1-866-641-0419 or by writing at the following address:

**Adapt Solutions/Warranty Department
145 Damase-Breton
St-Lambert-de-Lauzon, Québec
G0S 2W0 Canada**

WARRANTY REGISTRATION

Adapt Solutions #1 Ltd located in St-Lambert-de-Lauzon, Québec, Canada, warrants this product to be free from any manufacturer's defect and workmanship for a period of three years from the date of purchase. This product has to be installed by an authorized dealer and has to be operated properly to be warranted. This warranty is limited to the original purchaser and does not cover defects on the vehicle. The warranty registration must be completed and returned within 30 days of purchase. This warranty does not cover abnormal wear or damages resulting from an accident, misuse, neglect of if the product has been altered in any way or aspect. Defective parts must be returned by an authorised dealer, who will have received a prior return authorization from Adapt Solutions #1 Ltd. All warranted parts shall be sent to Adapt Solutions #1 Ltd prepaid. Parts will be replaced or repaired at Adapt Solutions #1 Ltd.

Date installed: _____

Dealer: _____

Dealer signature: _____

TO KEEP FOR YOUR RECORDS

(Place serial number here)

WE ARE THE SOLUTION

Please contact your dealer for immediate assistance

Dealer:

TF 1.866.641.0419 | tech@adaptsolutions.com



145 Damase-Breton | Saint-Lambert, Qc G0S 2W0

TF 1.866.641.0419 | **Fax** 418.889.9838 | **P** 418.889.0419

www.adaptsolutions.com

