



USER MANUAL



WE HAVE YOU COVERED

Our passionate team brings you innovative, quality mobility solutions. Our products are reliable, sleek, and designed with simplicity so that you can rest assured that they work when you need them.

Thanks to our out of the box thinking, we are vigilant to meet your mobility needs with our simple and practical solutions. We recognize the important role that safety plays in consumer choices and are proud to offer a line of products that maintain the safety features found in your vehicle. Backed by third party verification for compliance to applicable Federal Motor Vehicle Safety Standards and Canadian Motor Vehicle Safety Standards, our products are engineered for safety and durability.

Our commitment to quality and safety along with our genuine desire to meet your mobility needs makes for a product we are sure you will be proud to own.





XL-BASE

USER GUIDE

To get into your vehicle

The **XL-BASE** is a pivoting and lowering seat base that uses your vehicle's original seat which affords comfort and freedom. The user-friendly functions of the **XL-BASE** make it a breeze to operate.

Follow the next few steps:

- After opening the door, retrieve the hand held pendant.
- Press and hold the **OUT** button until the **XL-BASE** has reached the desired height for your ideal transfer.
- Once you are comfortably positioned on the **XL-BASE**, press and hold the **IN** button until the **XL-BASE** is completely locked in place. At any time you may release the **IN** button and the **XL-BASE** will stop to allow you to reposition. Once you are ready to continue, simply press and hold the **IN** button to resume your ascent into the vehicle.
- Close your vehicle door and replace the hand held pendant.



USER GUIDE (CONTINUED)

To get out of your vehicle

Follow the next few steps:

- Start by removing your seat belt and opening the door.
- Locate the hand held pendant.
- Press and hold the **OUT** button until the **XL-BASE** has reached the desired height for your ideal transfer. At any time you may release the **OUT** button and the **XL-BASE** will stop to allow you to reposition. Once you are ready to continue, simply press and hold the **OUT** button to resume you descent out of the vehicle.
- Once you have transferred out of the seat, press and hold the **IN** button until the **XL-BASE** has completely entered and locked in place, inside the vehicle.
- Replace the hand held pendant and close your vehicle door.

WARNING

NEVER TRY TO OPERATE THE **XL-BASE** WHILE THE VEHICLE DOOR IS CLOSED OR IF THE VEHICLE IS IN MOTION.

⚠ WARNING: This product can expose you to chemicals including benzene, which is known to the State of California to cause cancer and birth defects or other reproductive harm. For more information, go to www.P65Warnings.ca.gov.





MANUAL BACKUP SYSTEM



Your new **XL-BASE** is equipped with a manual backup system. If you were to encounter a power failure or other functional problem, this backup system will allow your seat to be returned to the inside of your vehicle.

PLEASE NOTE: This feature cannot be used with a passenger seated on the unit.

Please look through the follow pages for the steps to follow and details on how to use this feature. Please contact your dealer for full diagnostic and needed repairs.

MANUAL BACKUP SYSTEM

- If your **XL-BASE** will not run and the seat is **OUT** of the vehicle, follow the next few steps to manually raise the seat and pivot it back in.

Please note: The next few steps will require physical effort.

- Remove the vinyl covered trim that is located on the right-hand side of the **XL-BASE**. Under this trim, you will find two (2) 'Allen' keys (a red one and a yellow one).



MANUAL BACKUP SYSTEM (CONTINUED)

- Locate the yellow drive nut on the exterior of the slide (see picture for details).
- Using the biggest 'Allen' key (yellow 'Allen' key), unbolt the drive nut.
- Once you have removed this, you will be able to push the seat base up to its pivot point.



MANUAL BACKUP SYSTEM (CONTINUED)


- If the **pivot** function is defective, follow the next few steps.

Please note: The next few steps will require physical effort.

- Remove the vinyl covered trim that is located on the right-hand side of the **XL-BASE**. Under this trim you will find two (2) 'Allen' keys (a red one and a yellow one).
- Locate the red 'Allen' key. Use it to unbolt the red bolt holding the end of the actuator. Start by unbolting the actuator from the left end (motor end). Manually pivot the seat in about half way and then unbolt the other end of the actuator using the same 'Allen' key.



MANUAL BACKUP SYSTEM (CONTINUED)



***NOTE:** Be careful not to disconnect or constrict the wires when manually pivoting the seat. Once the problem is resolved, do not forget to replace the 'Allen' keys back in place, under the vinyl covered trim. Please contact your dealer for full diagnostic and needed repairs.

3 YEAR LIMITED WARRANTY XL-BASE

Adapt Solutions # 1 Ltd. (Adapt Solutions), warrants to the original purchaser of an **XL-BASE** that the equipment is free from defects in material and workmanship for a period of three years from date of purchase.

During the first year of the warranty, Adapt Solutions will supply the replacement parts as well as a pre-set monetary amount (determined by Adapt Solutions) for the repair if a defect in material or workmanship is discovered. After the initial year of this warranty, only parts and components are covered. This warranty does not cover labour and other services after the initial year. Freight and other related repair charges will be the responsibility of the original purchaser.

The only remedy for a defect in one of Adapt Solutions products (**XL-BASE; XL-SEAT; HI-LIFT; SPEEDY-LIFT; XL-BOARD; POWER-PULL**) shall be the repair or the replacement, at the discretion of Adapt Solutions, of the defective part or component. If repair or replacement is not commercially practical or cannot be timely made, Adapt Solutions may decide to refund the purchase price of the equipment instead of repairing or replacing the original equipment.

In no event shall Adapt Solutions be responsible for indirect, incidental or consequential damages, whether such damages arise from claims based on contract, warranty, tort (including negligence), strict liability or product liability.

All implied warranties, including any warranty of merchantability or fitness for a particular purpose, are limited in their duration to the length of the warranty stated above for the affected component.

This warranty is to the original purchaser only, and excludes product damage due to installation error, product misuse, product abuse, accidents, physical damage, damage in shipment, modifications not made by Adapt Solutions, or repairs undertaken by anyone other than authorized distributors.

TO OBTAIN WARRANTY PARTS OR REIMBURSEMENT

You must follow these procedures:

1. Obtain warranty authorization by calling your local Adapt Solutions dealer or Adapt Solutions directly at 1-866-641-0419.
2. Return the faulty Adapt Solutions component/equipment, freight prepaid, to the address provided by your Adapt Solutions dealer or Adapt Solutions with proof of purchase indicating the date purchased.

Adapt Solutions will pay for shipping back to the purchaser within the continental United States and Canada if a defect in materiel or workmanship is discovered. Return freight and repair charges will be the responsibility of the purchaser if the problem is not covered by warranty.

This warranty gives you specific legal rights, and you may also have other rights that vary from state to state/ province to province. Adapt Solutions specifically does not authorize any person to extend the time or scope of this warranty.

For further information regarding this limited warranty, please contact us by calling 1-866-641-0419 or by writing at the following address:

**Adapt Solutions / Warranty Department
145 Damase-Breton
St-Lambert-de-Lauzon, Québec
G0S 2W0 Canada**

WARRANTY REGISTRATION

Adapt Solutions #1 Ltd located in St-Lambert-de-Lauzon, Québec, Canada, warrants this product to be free from any manufacturer's defect and workmanship for a period of three years from the date of purchase. This product has to be installed by an authorized dealer and has to be operated properly to be warranted. This warranty is limited to the original purchaser and does not cover defects on the vehicle. The warranty registration must be completed and returned within 30 days of purchase. This warranty does not cover abnormal wear or damages resulting from an accident, misuse, neglect or if the product has been altered in any way or aspect. Defective parts must be returned by an authorized dealer, who will have received a prior return authorization from Adapt Solutions #1 Ltd. All warranted parts shall be sent to Adapt Solutions #1 Ltd prepaid. Parts will be replaced or repaired at Adapt Solutions #1 Ltd.

Date installed: _____

Dealer: _____

Dealer signature: _____

TO KEEP FOR YOUR RECORDS

(Place serial number here)

WE ARE THE SOLUTION

Please contact your dealer for immediate assistance

Dealer:

TF 1.866.641.0419 | tech@adaptsolutions.com



145 Damase-Breton | Saint-Lambert, Qc G0S 2W0

TF 1.866.641.0419 | **Fax** 418.889.9838 | **P** 418.889.0419

www.adaptsolutions.com

