

How it works

The XL-BOARD is a stationary, foldable and retractable vehicle transfer board designed to bridge the gap when transferring from a wheelchair to a driver or front passenger seat of the Mazda 5®.

Features & Specifications

- Assists transfer into the vehicle
- Easy to use
- Takes up very little space
- Retractable support arm
- Maximum capacity: 300 lbs
- Transfer board dimensions: 14.5" × 13"
- Transfer board height: 21.5" from the ground
- Mazda 5 seat height: 23" from the ground

It's a breeze

No structural modification to the vehicle means a quick and easy installation. The XL-BOARD installs directly into the seat's original floor bolts. It is also easily removed to return the vehicle to its original condition without affecting the resale value.



Mazda and Mazda 5 are registered trademarks of Mazda Motor Corporation





XL-BOARD

Transfer board for Mazda 5®

Great alternative

The XL-BOARD is designed to fit into the Mazda 5® because it is one of the only vehicles of its kind to offer a sliding door. This makes it an ideal feature for wheelchair users looking for an alternative to a minivan.

Join forces

Combining the XL-BOARD with the SPEEDY-LIFT may be the solution to simplify your vehicle access. It is important to note that a power sliding door is not available from Mazda or Adapt Solutions.

Adaptive equipment

Our products qualify as adaptive equipment and are eligible for reimbursement under the terms of Mobility programs offered by various vehicle manufacturers.

About us



We take our customers' safety very seriously. We are dedicated to manufacturing quality mobility equipment. Our products undergo a strict certification process.

Certification of our products involves different tests such as: crash testing, cycle testing and weight capacity testing. We are proud to offer a line of products that exceed the industry's standards.



Dealer:



145 Damase-Breton | Saint-Lambert, Qc G0S 2W0

TF 1.866.641.0419 | Fax 418.889.9838 | P 418.889.0419

www.adaptsolutions.com



Due to the diversity of mobility solutions available to you, please contact your local mobility dealer for a personalized evaluation and overview of options tailored to your needs.